

KERN 323-02

KERN

Weight set, F1, 1 mg - 50 g, Knob, stainless steel polished, in wooden case



Shipping weight 0,78 kg

Optional cases

Case - Wood article number 313-020-100

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Test weight - nominal value	1 mg - 50 g
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	152×152×66 mm
Shape weight	Knob
Material	stainless steel polished
Number of weights	20
Case of weight set	in wooden case

Services (optional)

Article number for conformity assessment	952-401
Article number for DAkkS calibration	962-401
Article number for guaranteed plus-deviation	921-401

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	150×140×65 mm
Shipping method	Parcel service
Net weight approx.	0,80 kg
Gross weight approx.	0,80 kg