

KERN 323-026

KERN

Weight set, F1, 1 mg - 50 g, Knob, stainless steel polished, in aluminium case



Shipping weight 0,8 kg

Optional cases

Case - Aluminium article number 313-020-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 mg - 50 g
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	165×120×75 mm
Shape weight	Knob
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for conformity assessment	952-401
Article number for DAkkS calibration	962-401
Article number for guaranteed plus-deviation	921-401

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Shipping method	Parcel service
Net weight approx.	0,65 kg
Gross weight approx.	0,80 kg