

# KERN 323-03

# KERN

Weight set, F1, 1 mg - 100 g, Knob, stainless steel polished, in wooden case



Shipping weight 0,88 kg

## Optional cases

Case - Wood article number 313-030-100

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

## Measuring System

Test weight - nominal value	1 mg - 100 g
OIML-compliant	✓
Weight OIML class	F1

## Construction

Dimension (W×D×H)	152×152×66 mm
Shape weight	Knob
Material	stainless steel polished
Number of weights	21
Case of weight set	in wooden case

## Services (optional)

Article number for conformity assessment	952-402
Article number for DAkkS calibration	962-402
Article number for guaranteed plus-deviation	921-402

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg