

# KERN 323-046

Weight set, F1, 1 mg - 200 g, Knob, stainless steel polished, in aluminium case



Shipping weight 1,1 kg

## Optional cases

Case - Aluminium article number 313-040-600

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

## Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 mg - 200 g
OIML-compliant	✓
Weight OIML class	F1

## Construction

Dimension (W×D×H)	190×130×80 mm
Shape weight	Knob
Material	stainless steel polished
Case of weight set	in aluminium case

## Services (optional)

Article number for conformity assessment	952-403
Article number for DAKKS calibration	962-403
Article number for guaranteed plus-deviation	921-403

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,2 kg