

KERN 323-064



Weight set, F1, 1 mg - 1 kg, Knob, stainless steel polished, in plastic case



Net weight approx.	3,8 kg
Gross weight approx.	3,8 kg
Shipping weight	4 kg
Shipping volume weight	4 kg

Optional cases

Case - Plastic article number	313-082-400
-------------------------------	-------------

Product information

GTIN/EAN number	4045761080713
Trade number	84239090
Warranty Time	36 mon
REACH classification	no prohibited substances

Pictograms

STANDARD



OPTION



Date of creation: 23.10.2026

Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	F1
OIML-compliant	✓

Construction

Dimension (W×D×H)	350×320×180 mm
Shape weight	Knob
Material	stainless steel polished
Number of weights	25
Case of weight set	in plastic case
is weight set	✓

Services (optional)

Article number for DAkkS calibration	962-405
Article number for guaranteed plus-deviation	921-405

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Shipping method	Parcel service