

KERN 323-066

KERN

Weight set, F1, 1 mg - 1 kg, Knob, stainless steel polished, in aluminium case



Shipping weight 3,5 kg

Optional cases

Case - Aluminium article number 313-060-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 mg - 1 kg
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	230×140×110 mm
Shape weight	Knob
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for conformity assessment	952-405
Article number for DAkkS calibration	962-405
Article number for guaranteed plus-deviation	921-405

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,6 kg