

KERN 323-086

KERN

Weight set, F1, 1 mg - 5 kg, Knob, stainless steel polished, in aluminium case



Shipping weight 18 kg

Optional cases

Case - Aluminium article number 313-080-600

Pictograms

STANDARD



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML F1 |
| Product family | 323 |

Measuring System

| | |
|--------------------------------|-------------|
| Conformity assessment possible | ✓ |
| Test weight - nominal value | 1 mg - 5 kg |
| OIML-compliant | ✓ |
| Weight OIML class | F1 |

Construction

| | |
|--------------------|--------------------------|
| Dimension (W×D×H) | 460×260×200 mm |
| Shape weight | Knob |
| Material | stainless steel polished |
| Case of weight set | in aluminium case |

Services (optional)

| | |
|--|---------|
| Article number for conformity assessment | 952-407 |
| Article number for DAKKS calibration | 962-407 |
| Article number for guaranteed plus-deviation | 921-407 |

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 460×260×200 mm |
| Shipping method | Parcel service |
| Net weight approx. | 13 kg |
| Gross weight approx. | 18 kg |