

Product brochure

08.08.2026

# KERN 323-096

Weight set, F1, 1 mg - 10 kg, Knob, stainless steel polished, in aluminium case

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

☎ +49 7433 9933-0  
✉ [info@kern-sohn.com](mailto:info@kern-sohn.com)  
🌐 [www.kern-sohn.com](http://www.kern-sohn.com)



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

| Model   | Description                                          |
|---------|------------------------------------------------------|
| 921-408 | Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-408 |
| 970-112 | Warranty extension (+2 years) 970-112                |

SERVICES

| Model    | Description                                          |
|----------|------------------------------------------------------|
| 920-709  | Marking (Last digits of serial number) KERN 920-709  |
| 921-408  | Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-408 |
| 962-408  | DAkKS calibration certificate 962-408                |
| 962-408R | DAkKS Recalibration 962-408R                         |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days