

KERN 324-046

Weight set, F1, 1 g - 200 g, Knob, stainless steel polished, in aluminium case



Shipping weight 1,1 kg

Optional cases

Case - Aluminium article number 314-040-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	324

Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 g - 200 g
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	190×130×80 mm
Shape weight	Knob
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for conformity assessment	952-417
Article number for DAkkS calibration	962-417
Article number for guaranteed plus-deviation	921-417

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	1,2 kg