

Product brochure

08.08.2026

KERN 325-062

Weight set, F1, 1 mg - 1 kg, Eco Shape, stainless steel polished, in wooden case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
313-060-100	Wood weight case 313-060-100 (for nominal values 1 mg - 1 kg for classes E1+E2+F1 for design Cylindrical/compact)

SERVICES

Model	Description
920-706	Marking (Last digits of serial number) KERN 920-706
921-405	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-405
962-405	DAkKS calibration certificate 962-405
962-405R	DAkKS Recalibration 962-405R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days