

Product brochure

08.08.2026

# KERN 325-072

Weight set, F1, 1 mg - 2 kg, Eco Shape, stainless steel polished, in wooden case


ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

| Model       | Description                                                                                                           |
|-------------|-----------------------------------------------------------------------------------------------------------------------|
| 313-070-100 | Wood weight case 313-070-100 (for nominal values 1 mg - 2 kg   for classes E1+E2+F1   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                             |
|----------|---------------------------------------------------------|
| 920-707  | Marking (Last digits of serial number) KERN 920-707     |
| 921-406  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-406 |
| 962-406  | DAkkS calibration certificate 962-406                   |
| 962-406R | DAkkS Recalibration 962-406R                            |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days