

Product brochure

08.08.2026

KERN 325-096

Weight set, F1, 1 mg - 10 kg, Eco Shape, stainless steel polished, in aluminium case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-090-600	Aluminium weight case 313-090-600 (for nominal values 1 mg - 10 kg for classes E1 - M1 for design Cylindrical/compact)

SERVICES

Model	Description
920-709	Marking (Last digits of serial number) KERN 920-709
921-408	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-408
962-408	DAkkS calibration certificate 962-408
962-408R	DAkkS Recalibration 962-408R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days