

# KERN 326-01

# KERN

Individual weight, OIML Class F1, 1 g, Eco Shape, stainless steel polished



Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	347-030-400
Case - Aluminium article number	317-010-600
Case - Wood article number	317-010-100

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

## Measuring System

Test weight - Density (kg/m <sup>3</sup> )	7950
Test weight - Density uncertainty (kg/m <sup>3</sup> )	140
Test weight - Density [Min] (OIML)	3000
Test weight - nominal value	1 g
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	0,1 mg

## Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	6×6 mm
Material	stainless steel polished

## Services (optional)

Article number for DAkkS calibration	962-431
Article number for guaranteed plus-deviation	921-431

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	6×6×6 mm
Shipping method	Parcel service