

# KERN 326-022

# KERN

Weight set, F1, 1 g - 50 g, Eco Shape, stainless steel polished, in wooden case



## Optional cases

Case - Wood article number 313-020-100

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	326

## Measuring System

Test weight - nominal value	1 g - 50 g
OIML-compliant	✓
Weight OIML class	F1

## Construction

Dimension (W×D×H)	152×152×66 mm
Shape weight	Eco Shape
Material	stainless steel polished
Case of weight set	in wooden case

## Services (optional)

Article number for DAkkS calibration	962-415
Article number for guaranteed plus-deviation	921-415

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Shipping method	Parcel service
Net weight approx.	0,80 kg
Gross weight approx.	0,80 kg
Shipping weight	0,78 kg