

KERN 326-026

KERN

Weight set, F1, 1 g - 50 g, Eco Shape, stainless steel polished, in aluminium case



Optional cases

Case - Aluminium article number 314-020-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	1 g - 50 g
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	165×120×75 mm
Shape weight	Eco Shape
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	962-415
Article number for guaranteed plus-deviation	921-415

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg