

Product brochure

07.08.2026

KERN 326-02

Individual weight, OIML Class F1, 2 g, Eco Shape, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
317-020-100	Wood weight case 317-020-100 (for nominal values 2 g for classes E1+E2+F1 for design Cylindrical/compact)
317-020-600	Aluminium weight case 317-020-600 (for nominal values 2 g for classes E1 - M3 for design Cylindrical/compact)
347-030-400	Plastic weight case 347-030-400 (for nominal values 1 g - 5 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-432	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-432
962-432	DAkks calibration certificate 962-432
962-432R	DAkks Recalibration 962-432R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days