

Product brochure

08.08.2026

KERN 326-034

Weight set, F1, 1 g - 100 g, Eco Shape, stainless steel polished, in plastic case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
3 13-052-400	Plastic carrying case 3 13-052-400 (for nominal values bis 500 g for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-711	Marking (Last digits of serial number) KERN 920-711
921-416	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-416
962-416	DAkkS calibration certificate 962-416
962-416R	DAkkS Recalibration 962-416R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days