

KERN 326-042

KERN

Weight set, F1, 1 g - 200 g, Eco Shape, stainless steel polished (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	326

Measuring System

OIML-compliant	✓
Weight OIML class	F1
Test weight - nominal value	1 g - 200 g

Construction

Shape weight	Eco Shape
Dimension (W×D×H)	152×152×66 mm
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Net weight	0,78 kg
Shipping method	Parcel service
Net weight approx.	0,80 kg
Gross weight approx.	0,80 kg
Shipping weight	0,78 kg

Services (optional)

Article number for DAkkS calibration	962-417
Article number for guaranteed plus-deviation	921-417