

Product brochure

08.08.2026

KERN 326-052

Weight set, F1, 1 g - 500 g, Eco Shape, stainless steel polished, in wooden case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-050-100	Wood weight case 313-050-100 (for nominal values 1 mg - 500 g for classes E1+E2+F1 for design Cylindrical/compact)

SERVICES

Model	Description
920-713	Marking (Last digits of serial number) KERN 920-713
921-418	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-418
962-418	DAkkS calibration certificate 962-418
962-418R	DAkkS Recalibration 962-418R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days