

Product brochure

08.08.2026

# KERN 326-054

Weight set, F1, 1 g - 500 g, Eco Shape, stainless steel polished, in plastic case

ENGLISH




## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com



PICTOGRAMS



ACCESSORIES

| Model       | Description                                                                                                             |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| 313-052-400 | Plastic carrying case 313-052-400 (for nominal values bis 500 g   for classes E1 - M3   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                             |
|----------|---------------------------------------------------------|
| 920-713  | Marking (Last digits of serial number) KERN 920-713     |
| 921-418  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-418 |
| 962-418  | DAkKS calibration certificate 962-418                   |
| 962-418R | DAkKS Recalibration 962-418R                            |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days