

Product brochure

08.08.2026

KERN 326-064

Weight set, F1, 1 g - 1 kg, Eco Shape, stainless steel polished, in plastic case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-082-400	Plastic carrying case 313-082-400 (for nominal values bis 5 kg for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-714	Marking (Last digits of serial number) KERN 920-714
921-419	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-419
962-419	DAkkS calibration certificate 962-419
962-419R	DAkkS Recalibration 962-419R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days