

Product brochure

07.08.2026

KERN 326-06

Individual weight, OIML Class F1, 50 g, Eco Shape, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
317-060-100	Wood weight case 317-060-100 (for nominal values 50 g for classes E1+E2+F1 for design Cylindrical/compact)
317-060-600	Aluminium weight case 317-060-600 (for nominal values 50 g for classes E1 - M3 for design Cylindrical/compact)
347-070-400	Plastic weight case 347-070-400 (for nominal values 50 g - 100 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-436	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-436
962-436	DAkkS calibration certificate 962-436
962-436R	DAkkS Recalibration 962-436R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days