

Product brochure

08.08.2026

KERN 326-074

Weight set, F1, 1 g - 2 kg, Eco Shape, stainless steel polished, in plastic case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-082-400	Plastic carrying case 313-082-400 (for nominal values bis 5 kg for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-715	Marking (Last digits of serial number) KERN 920-715
921-420	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-420
962-420	DAkkS calibration certificate 962-420
962-420R	DAkkS Recalibration 962-420R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days