

Product brochure

08.08.2026

KERN 326-14

Individual weight, OIML Class F1, 10 kg, Eco Shape, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
317-140-100	Wood weight case 317-140-100 (for nominal values 10 kg for classes E1+E2+F1 for design Cylindrical/compact)
317-140-600	Aluminium weight case 317-140-600 (for nominal values 10 kg for classes E1 - M3 for design Cylindrical/compact)
347-140-400	Plastic weight case 347-140-400 (for nominal values 10 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-444	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-444
962-444	DAkkS calibration certificate 962-444
962-444R	DAkkS Recalibration 962-444R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days