

Product brochure

08.08.2026

KERN 326-36

Individual weight, OIML Class F1, 5 kg, Block, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
346-060-600	Aluminium weight case 346-060-600 (for nominal values 5 kg Block for classes F1 - M3 for design Block)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-443	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-443
962-443	DAkkS calibration certificate 962-443
962-443R	DAkkS Recalibration 962-443R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days