

Product brochure

08.08.2026

# KERN 326-37

Individual weight, OIML Class F1, 10 kg, Block, stainless steel polished

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

## PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model       | Description                                                                                                 |
|-------------|-------------------------------------------------------------------------------------------------------------|
| 346-070-600 | Aluminium weight case 346-070-600 (for nominal values 10 kg Block   for classes F1 - M3   for design Block) |

## SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 921-444  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-444             |
| 952-444R | (Re-)Calibration 952-444R                                           |
| 962-444  | DAkks calibration certificate 962-444                               |
| 962-444R | DAkks Recalibration 962-444R                                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days