

KERN 326-37

KERN

Individual weight, OIML Class F1, 10 kg, Block, stainless steel polished



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Conformity assessment possible	✓
Test weight - volume uncertainty (cm ³)	15,625
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Test weight - nominal value	10 kg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	50 mg

Construction

Shape weight	Block
Dimensions weight (W×D×H)	193×96×117 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-444
Article number for DAkkS calibration	962-444

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	193×96×117 mm
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg

Optional cases

Case - Plastic article number	nicht verfügbar
Case - Aluminium article number	346-070-600
Case - Wood article number	nicht verfügbar

Pictograms

STANDARD

