

Product brochure

08.08.2026

KERN 326-39

Individual weight, OIML Class F1, 50 kg, Block, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
346-090-600	Aluminium weight case 346-090-600 (for nominal values 50 kg Block for classes F1 - M3 for design Block)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-446	DAkkS calibration certificate 962-446
962-446R	DAkkS Recalibration 962-446R

KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days