

Product brochure

07.08.2026

KERN 327-03

Individual weight, OIML Class F1, 5 g, Knob, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
315-250	Tweezers 315-250 (stainless steel 150 mm)
317-030-100	Wood weight case 317-030-100 (for nominal values 5 g for classes E1+E2+F1 for design Cylindrical/compact)
317-030-600	Aluminium weight case 317-030-600 (for nominal values 5 g for classes E1 - M3 for design Cylindrical/compact)
347-030-400	Plastic weight case 347-030-400 (for nominal values 1 g - 5 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-433	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-433
962-433	DAkkS calibration certificate 962-433
962-433R	DAkkS Recalibration 962-433R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days