

# KERN 327-03

# KERN

Individual weight, OIML Class F1, 5 g, Knob, stainless steel polished (OIML)



Article number for DAkkS calibration 962-433

Article number for guaranteed plus-deviation 921-433

## Packing & Shipping

Delivery time 1 d  
Dimensions packaging (W×D×H) 8×8×16 mm  
Net weight 0,005 kg  
Shipping method Parcel service  
Net weight approx. 0,05 kg  
Gross weight approx. 0,05 kg  
Shipping weight 0,1 kg

## Optional cases

Case - Plastic article number 347-030-400  
Case - Aluminium article number 317-030-600  
Case - Wood article number 317-030-100

## Pictograms

## STANDARD



## Category

Brand KERN  
Product category Individual weight  
Product group OIML F1  
Product family 327

## Measuring System

Conformity assessment possible ✓  
Test weight - volume (cm<sup>3</sup>) 0,625  
Test weight - volume uncertainty (cm<sup>3</sup>) 0,0078  
Test weight - Density (kg/m<sup>3</sup>) 7950  
Test weight - Density uncertainty (kg/m<sup>3</sup>) 140  
Test weight - Density [Min] (OIML) 5300  
Test weight - Density [Max] (OIML) 16000  
Test weight - nominal value 5 g  
OIML-compliant ✓  
Weight OIML class F1  
Permissible tolerance according to OIML (+/-) 0,16 mg

## Construction

Shape weight Knob  
Dimensions (Ø×H) 8×16,5 mm  
Material stainless steel polished (OIML)

## Services (optional)

Article number for conformity assessment 952-433