

Product brochure

07.08.2026

KERN 327-141

Individual weight, OIML Class F1, 10 kg, Check weight, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
337-141-100	Wood weight case 337-141-100 (for nominal values 10 kg for classes F1 for design Check weight)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-444	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-444
962-444	DAkKS calibration certificate 962-444
962-444R	DAkKS Recalibration 962-444R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days