

Product brochure

07.08.2026

# KERN 327-151

Individual weight, OIML Class F1, 20 kg, Check weight, stainless steel polished

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option





ACCESSORIES

| Model       | Description                                                                                        |
|-------------|----------------------------------------------------------------------------------------------------|
| 337-151-100 | Wood weight case 337-151-100 (for nominal values 20 kg   for classes F1   for design Check weight) |

SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 962-445  | DAkkS calibration certificate 962-445                               |
| 962-445R | DAkkS Recalibration 962-445R                                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days