

Product brochure

07.08.2026

KERN 327-161

Individual weight, OIML Class F1, 50 kg, Check weight, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 +49 7433 9933-0
 info@kern-sohn.com
 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
337-161-100	Wood weight case 337-161-100 (for nominal values 50 kg for classes F1 for design Check weight)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-446	DAkkS calibration certificate 962-446
962-446R	DAkkS Recalibration 962-446R

KERN Pictogram Overview

	Pallet shipment: The duration of internal product provision are two days
	Accredited calibration (DKD): The duration of the accredited calibration is 3 working days