

# KERN 328-222-410

**KERN**

Plastic box for weight sets (Cylindrical/compact, 1 mg - 5 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

## Measuring System

Case suitable for OIML classes	E1 - M2
Case suitable for nominal weight values	1 mg - 5 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	130×60×22 mm
Material	plastic
Case - Smallest possible weight	1 mg
Case - Highest possible weight	5 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	130×60×22 mm
Net weight	0,09 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

## Pictograms