

Product brochure

08.08.2026

KERN 328-222-420

Plastic box for weight sets (Cylindrical/compact, 1 mg - 5 g)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
318-270	Brush 318-270

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day