

Product brochure

08.08.2026

KERN 328-226

Set of milligram weights, F1, 1 mg - 500 mg, Platelet, stainless steel polished, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
315-247	Tweezers 315-247 (stainless steel, plastic 136 mm)

SERVICES

Model	Description
920-701	Marking (Last digits of serial number) KERN 920-701
921-450	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-450
962-450	DAkkS calibration certificate 962-450
962-450R	DAkkS Recalibration 962-450R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days