

KERN 328-226

KERN

Set of milligram weights, F1, 1 mg - 500 mg, Platelet, stainless steel polished, in aluminium case



Optional cases

Case - Aluminium article number 313-010-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	328

Measuring System

Test weight - nominal value	1 mg - 500 mg
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	165×158×37 mm
Shape weight	Platelet
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	962-450
Article number for guaranteed plus-deviation	921-450

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×158×37 mm
Shipping method	Parcel service
Net weight approx.	0,40 kg
Gross weight approx.	0,40 kg
Shipping weight	0,39 kg