

Product brochure

07.08.2026

KERN 333-034

Weight set, F2, 1 mg - 100 g, Knob, stainless steel fine turned, in plastic case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-052-400	Plastic carrying case 313-052-400 (for nominal values bis 500 g for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-703	Marking (Last digits of serial number) KERN 920-703
921-402	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-402
962-402	DAkKS calibration certificate 962-402
962-402R	DAkKS Recalibration 962-402R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days