

Product brochure

07.08.2026

KERN 333-04

Weight set, F2, 1 mg - 200 g, Knob, stainless steel fine turned, in wooden case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
353-420-210	Wood weight case 353-420-210

SERVICES

Model	Description
920-704	Marking (Last digits of serial number) KERN 920-704
921-403	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-403
962-403	DAkKS calibration certificate 962-403
962-403R	DAkKS Recalibration 962-403R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days