

Product brochure

07.08.2026

KERN 333-054

Weight set, F2, 1 mg - 500 g, Knob, stainless steel fine turned, in plastic case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
313-052-400	Plastic carrying case 313-052-400 (for nominal values bis 500 g for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-705	Marking (Last digits of serial number) KERN 920-705
921-404	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-404
962-404	DAkKS calibration certificate 962-404
962-404R	DAkKS Recalibration 962-404R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days