

Product brochure

08.08.2026

KERN 333-066

Weight set, F2, 1 mg - 1 kg, Knob, stainless steel fine turned, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany


☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS



ACCESSORIES

Model	Description
313-060-600	Aluminium weight case 313-060-600 (for nominal values 1 mg - 1 kg for classes E1 - M1 for design Cylindrical/compact)
313-082-400	Plastic carrying case 313-082-400 (for nominal values bis 5 kg for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-706	Marking (Last digits of serial number) KERN 920-706
921-405	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-405
962-405	DAkKS calibration certificate 962-405
962-405R	DAkKS Recalibration 962-405R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days