

Product brochure

07.08.2026

KERN 333-09

Weight set, F2, 1 mg - 10 kg, Knob, stainless steel fine turned, in wooden case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
353-490-200	Wood weight case 353-490-200 (for nominal values 1 mg - 10 kg for classes F2 + M1 + M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-709	Marking (Last digits of serial number) KERN 920-709
921-408	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-408
962-408	DAkKS calibration certificate 962-408
962-408R	DAkKS Recalibration 962-408R

KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days