

Product brochure

08.08.2026

KERN 334-036

Weight set, F2, 1 g - 100 g, Knob, stainless steel fine turned, in aluminium case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
314-030-600	Aluminium weight case 314-030-600 (for nominal values 1 g - 100 g for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-711	Marking (Last digits of serial number) KERN 920-711
921-416	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-416
962-416	DAkKS calibration certificate 962-416
962-416R	DAkKS Recalibration 962-416R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days