

Product brochure

07.08.2026

KERN 334-04

Weight set, F2, 1 g - 200 g, Knob, stainless steel fine turned, in wooden case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
353-440-200	Wood weight case 353-440-200 (for nominal values 1 mg - 200 g for classes F2 + M1 + M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-712	Marking (Last digits of serial number) KERN 920-712
921-417	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-417
962-417	DAkKS calibration certificate 962-417
962-417R	DAkKS Recalibration 962-417R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days