

# KERN 334-056

**KERN**

Weight set, F2, 1 g - 500 g, Knob, stainless steelfine turned, in aluminium case



## Optional cases

Case - Aluminium article number 314-050-600

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

## Measuring System

Test weight - nominal value	1 g - 500 g
OIML-compliant	✓
Weight OIML class	F2

## Construction

Dimension (W×D×H)	210×140×98 mm
Shape weight	Knob
Material	stainless steelfine turned
Case of weight set	in aluminium case

## Services (optional)

Article number for DAkkS calibration	962-418
Article number for guaranteed plus-deviation	921-418

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	210×140×98 mm
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg