

Product brochure

07.08.2026

KERN 334-06

Weight set, F2, 1 g - 1 kg, Knob, stainless steel fine turned, in wooden case



ENGLISH

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option





SERVICES

Model	Description
920-714	Marking (Last digits of serial number) KERN 920-714
921-419	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-419
962-419	DAkKS calibration certificate 962-419
962-419R	DAkKS Recalibration 962-419R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days