

Product brochure

08.08.2026

KERN 334-066

Weight set, F2, 1 g - 1 kg, Knob, stainless steel fine turned, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
314-060-600	Aluminium weight case 314-060-600 (for nominal values 1 g - 1 kg for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-714	Marking (Last digits of serial number) KERN 920-714
921-419	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-419
962-419	DAkkS calibration certificate 962-419
962-419R	DAkkS Recalibration 962-419R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days