

KERN 335-070-200

KERN

F2 + M3 Wooden box for individual weight sets 1 mg - 2 kg



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for nominal weight values	1 mg - 2 kg
Case suitable for OIML classes	F2 + M3
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	305×180×140 mm
Material	beech
Case - Highest possible weight	2 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	310×180×140 mm
Shipping method	Parcel service
Net weight approx.	3,0 kg
Gross weight approx.	3,0 kg
Shipping weight	2,9 kg

Product information

GTIN/EAN number	4045761006140
Trade number	84239090
REACH classification	no prohibited substances