

# KERN 337-030-200

**KERN**

Wooden box for individual weights (Cylindrical/compact, 5 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	5 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	60×60×45 mm
Material	beech
Case - Smallest possible weight	5 g
Case - Highest possible weight	5 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	60×60×45 mm
Net weight	0,12 kg
Shipping method	Parcel service
Net weight approx.	0,15 kg
Gross weight approx.	0,15 kg
Shipping weight	0,12 kg

## Pictograms

STANDARD

