

Product brochure

07.08.2026

# KERN 337-04

Individual weight, OIML Class F2, 10 g, Knob, stainless steel fine turned

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

## PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model       | Description                                                                                                             |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| 315-242     | Tweezers 315-242 (polyamide   100 mm)                                                                                   |
| 315-243     | Tweezers 315-243 (stainless steel   105 mm)                                                                             |
| 315-250     | Tweezers 315-250 (stainless steel   150 mm)                                                                             |
| 317-040-600 | Aluminium weight case 317-040-600 (for nominal values 10 g   for classes E1 - M3   for design Cylindrical/compact)      |
| 337-040-200 | Wood weight case 337-040-200 (for nominal values 10 g   for classes F2 + M1   for design Cylindrical/compact)           |
| 347-050-400 | Plastic weight case 347-050-400 (for nominal values 10 g - 20 g   for classes F1 - M3   for design Cylindrical/compact) |

## SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 921-434  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-434             |
| 962-434  | DAkKS calibration certificate 962-434                               |
| 962-434R | DAkKS Recalibration 962-434R                                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days